Spectrum IQ Infusion System Integrating the Healthcare Enterprise (IHE) Integration Statement

The IHE Integration Statement below describes the intended conformance of the **Spectrum IQ** Infusion System with the IHE Technical Framework.

Vendor	Product Name		Date
Baxter	Spectrum IQ Gateway		January 2020
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Profiles Implemented	Actor Implemented		Options Implemented
Device Enterprise Communications (DEC)	Device Observation Reporter (DOR)		PCD-01 (Communicated PCD Data)
Infusion Pump Event Communications (IPEC)	IPEC - DOR (Device Observation Reporter)		PCD-10 (Communicate Infusion Pump Event)
Point-of-Care Infusion Verification (PIV)	Infusion Order Consumer		PCD-03 (Communicate Infusion Order)
Alert Communication Management (ACM)	Alert Reporter (AR)		PCD-04 (Report Alert)
Consistent Time (C1)	Time Client		ITI TF-1
Medical Equipment Management Location Services (MEMLS)	Location Observation Consumer (LOC)		PCD-16 (Report Location Observation)
Links to general information on IHE			
In North America: www.ihe.net		Profiles: www.ihe.net/Patient_Care_Devices/	

Rx Only. For safe and proper use of this device, refer to the appropriate operator's manuals.

www.baxter.com www.SpectrumIQ.com

Baxter Healthcare Corporation
One Baxter Parkway / Deerfield IL 60015
Baxter and Spectrum IQ are trademarks of Baxter International Inc.

