

Baxter

NovumIQ Syringe

INFUSION SYSTEM

COMPATIBLE SYRINGES LIST

COMPATIBLE WITH THE
NOVUM IQ SYRINGE PUMP



CODE	PRODUCT NAME	MIN RATE (ML/HR)	MAX RATE (ML/HR)
BD LUER-LOK			
309628	1 mL BD LUER-LOK SYRINGE STERILE, SINGLE USE POLYCARBONATE	0.01	30
309657	3 mL BD LUER-LOK SYRINGE STERILE, SINGLE USE	0.01	100
309646	5 mL BD LUER-LOK SYRINGE STERILE, SINGLE USE	0.03	200
309604	10 mL BD LUER-LOK SYRINGE STERILE, SINGLE USE	0.05	300
302830	20 mL BD LUER-LOK SYRINGE STERILE, SINGLE USE	0.1	500
302832	30 mL BD LUER-LOK SYRINGE STERILE, SINGLE USE	0.1	650
309653	50 mL BD LUER-LOK SYRINGE STERILE, SINGLE USE	0.1	1200
BD ENTERAL			
305853	3 mL SYRINGE ENTERAL 3 mL	0.01	100
305855	5 mL BD 5 mL ORAL ENTERAL SYRINGE	0.03	200
305857	10 mL SYRINGE ENTERAL 10 mL	0.05	300
305859	20 mL BD 20 mL ORAL ENTERAL SYRINGE	0.1	500
305863	50 mL SYRINGE ENTERAL 50 mL	0.1	1200

CODE	PRODUCT NAME	MIN RATE (ML/HR)	MAX RATE (ML/HR)
MONOJECT			
1180300777	SYRINGE 3 cc LUER LOCK TIP, WIDE FLANGE	0.01	105
1180600777 ¹	SYRINGE 6 CC LUER LOCK TIP, WIDE FLANGE	0.03	213
8881516937	SYRINGE 6 cc LUER LOCK TIP, STD FLANGE	0.03	213
1181200777 ¹	SYRINGE 12 cc LUER LOCK TIP, WIDE FLANGE	0.05	325
8881512878	SYRINGE 12 cc LUER LOCK TIP, STD FLANGE	0.05	325
1182000777, 8881520657	SYRINGE 20 cc LUER LOCK TIP	0.1	536
1183500777, 8881535762	SYRINGE 35 cc LUER LOCK TIP	0.1	735
1186000777	SYRINGE 60 cc LUER LOCK TIP, WIDE FLANGE	0.2	944
8881560125 ¹	SYRINGE 60 cc LUER LOCK TIP, STD FLANGE	0.2	944

BAXTER, ENFIT			
H93899106	6 mL ENTERAL SYRINGE WITH ENFIT CONNECTOR 6 mL, CLEAR WITH ORANGE PLUNGER	0.03	213
H93899112	12 mL ENTERAL SYRINGE WITH ENFIT CONNECTOR 12 mL, CLEAR WITH ORANGE PLUNGER	0.05	325
H93899120	20 mL ENTERAL SYRINGE WITH ENFIT CONNECTOR 20 mL, CLEAR WITH ORANGE PLUNGER	0.1	536
H93899160	60 mL ENTERAL SYRINGE WITH ENFIT CONNECTOR 60 mL, CLEAR WITH ORANGE PLUNGER	0.1	944

NEOMED, ENFIT			
NM-S6NC, PNM-S6NC	6 ML SYRINGE ORAL 6 ML WITH ENFIT	0.03	213
NM-S12NC, PNM-S12NC	12 mL SYRINGE ORAL 12 mL WITH ENFIT	0.05	325
NM-S20NC, PNM-S20NC	20 mL SYRINGE ORAL 20 mL WITH ENFIT	0.1	536
NM-S35NC, PNM-S35NC	35 mL SYRINGE ORAL 35 mL WITH ENFIT	0.1	735
NM-S60NC, PNM-S60NC	60 mL SYRINGE ORAL 60 mL NEOCONNECT WITH ENFIT	0.1	944

NOTE: ¹ Because of variations among models, these Monoject Syringes may in some cases not fully empty. The Syringe Empty Alarm occurs with some fluid left in the syringe.

This compatible syringe list is subject to change over time.

Please visit <https://service.baxter.com/documents> for the most current version.