

COMPATIBLE ADMINISTRATION SETS

COMPATIBLE WITH
 THE **NOVUM IQ**
 LARGE VOLUME PUMP



Compatible administration sets are designated as either PRIMARY compatible sets or OTHER compatible sets.

Volumetric accuracy for PRIMARY compatible sets is as stated below:

RATE	ACCURACY
0.5-0.9 mL/hr	±10%
01-1200 mL/hr	±5%

Volumetric accuracy for OTHER compatible sets is as stated below:

RATE	ACCURACY
0.5-1200 mL/hr	±10%

CODE BRIEF DESCRIPTION

BASIC SOLUTION SETS, 10 DROPS/ML

1C8109	SOLUTION SET, MALE LUER LOCK ADAPTER
1C8160	SOLUTION SET, MALE LUER LOCK ADAPTER
1C8296	SOLUTION SET, MALE LUER LOCK ADAPTER
1C8507	SOLUTION SET WITH DUO-VENT SPIKE, MALE LUER LOCK ADAPTER
1C8581	INTERLINK SYSTEM, Y-TYPE SOLUTION SET, 2 INJECTION SITES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C6401	INTERLINK SYSTEM, SOLUTION SET, INJECTION SITE, LEVER LOCK CANNULA, LUER LOCK ADAPTER
2C6419	INTERLINK SYSTEM, SOLUTION SET WITH DUO-VENT SPIKE, INJECTION SITE, MALE LUER LOCK ADAPTER

Rx only. For safe and proper use of the products mentioned herein, please refer to the appropriate Operator's Manual or Instructions for Use.

CODE	BRIEF DESCRIPTION
------	-------------------

BASIC SOLUTION SETS, 10 DROPS/ML (CONT.)

2C6425	INTERLINK SYSTEM, SOLUTION SET, 2 INJECTION SITES, MALE LUER LOCK ADAPTER
2C8401	CLEARLINK SYSTEM, SOLUTION SET, LUER ACTIVATED VALVE, MALE LUER LOCK
2C8419	CLEARLINK SYSTEM, SOLUTION SET WITH DUO-VENT SPIKE, LUER ACTIVATED VALVE, MALE LUER LOCK ADAPTER
2C8425	CLEARLINK SYSTEM, SOLUTION SET, 2 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER
ACT5435	SOLUTION SET WITH DUO-VENT SPIKE, INJECTION SITE, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
ACT5439	SOLUTION SET WITH DUO-VENT SPIKE, 2 INJECTION SITES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
ACT8419	CLEARLINK SYSTEM, SOLUTION SET WITH DUO-VENT SPIKE, LUER ACTIVATED VALVE, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR

BASIC SOLUTION SETS, 60 DROPS/ML

2C6402 ¹	INTERLINK SYSTEM, SOLUTION SET, INJECTION SITE, LEVER LOCK CANNULA MALE LUER LOCK ADAPTER
2C6424 ¹	INTERLINK SYSTEM, SOLUTION SET, 2 INJECTION SITES, MALE LUER LOCK ADAPTER
2C7554 ¹	SOLUTION SET, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR. FOR EPIDURAL USE
2C8402 ¹	CLEARLINK SYSTEM, SOLUTION SET, LUER ACTIVATED VALVE, MALE LUER LOCK ADAPTER
2C8554 ¹	SOLUTION SET, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
ACT5432 ¹	SOLUTION SET, INJECTION SITE, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR

BLOOD/SOLUTION SETS

1C8372	INTERLINK SYSTEM, Y-TYPE BLOOD/SOLUTION SET WITH STANDARD BLOOD FILTER, 170 TO 260 MICRON FILTER, 2 INJECTION SITES, MALE LUER LOCK ADAPTER
1C8732	Y-TYPE BLOOD/SOLUTION SET WITH STANDARD BLOOD FILTER, 170 TO 260 MICRON FILTER, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C6750H	INTERLINK SYSTEM, Y-TYPE BLOOD /SOLUTION SET, STANDARD BLOOD FILTER, 170 TO 260 MICRON FILTER, INJECTION SITE, MALE LUER LOCK ADAPTER
2C8750	CLEARLINK SYSTEM, Y-TYPE BLOOD/SOLUTION SET WITH STANDARD BLOOD FILTER, 170 TO 260 MICRON FILTER, LUER ACTIVATED VALVE, MALE LUER LOCK ADAPTER
2C8751	CLEARLINK SYSTEM, STRAIGHT TYPE BLOOD/SOLUTION SET WITH STANDARD BLOOD FILTER, 170 TO 260 MICRON FILTER, LUER ACTIVATED VALVE, MALE LOCK ADAPTER WITH RETRACTABLE COLLAR

BURETROL SETS

2C7562	INTERLINK SYSTEM, BURETROL SOLUTION SET, 150 ML VALVELESS BURETTE, 4 INJECTION SITES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C7564	INTERLINK SYSTEM, BURETROL SOLUTION SET, 150 ML BURETTE, DRIP CHAMBER FILTER VALVE, 3 INJECTION SITES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C8862	CLEARLINK SYSTEM, BURETROL SOLUTION SET, 150 ML VALVELESS BURETTE, 4 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C8864	CLEARLINK SYSTEM, BURETROL SOLUTION SET, 150 ML BURETTE, DRIP CHAMBER FILTER VALVE, 3 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR

CODE	BRIEF DESCRIPTION
------	-------------------

CONTINU-FLO SETS, 10 DROPS/ML

1C8722	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, 3 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER. EXTENSION SET, 2 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER. EXTENSION SET, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR.
1C8727	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, 3 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C6519	INTERLINK SYSTEM, CONTINU-FLO SOLUTION SET, 2 INJECTION SITES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C6537	INTERLINK SYSTEM, CONTINU-FLO SOLUTION SET, 3 INJECTION SITES, MALE LUER LOCK ADAPTER
2C6541	INTERLINK SYSTEM, CONTINU-FLO SOLUTION SET WITH DUO-VENT SPIKE, 3 INJECTION SITES, MALE LUER LOCK ADAPTER
2C8515	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, LUER ACTIVATED VALVE, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C8519	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, 2 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER
2C8537	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, 3 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C8541	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET WITH DUO-VENT SPIKE, 3 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
ACT5543	CONTINU-FLO SOLUTION SET WITH DUO-VENT SPIKE, 2 INJECTION SITES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR

CONTINU-FLO SETS, 60 DROPS/ML

2C6546¹	INTERLINK SYSTEM, CONTINU-FLO SOLUTION SET, 3 INJECTION SITES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C8546¹	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, 3 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C8548¹	CLEARLINK SYSTEM, VENTED CONTINU-FLO SOLUTION SET, 2 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR

FILTRATION SETS

2H6480^{1,2}	INTERLINK SYSTEM, NON-DEHP SOLUTION SET WITH DUO-VENT SPIKE, 0.2 MICRON DOWNSTREAM FILTER, INJECTION SITE, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C6571¹	INTERLINK SYSTEM, CONTINU-FLO SOLUTION SET (10 DROPS/ML), 0.2 MICRON HIGH PRESSURE, EXTENDED LIFE FILTER, 2 INJECTION SITES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C6572¹	INTERLINK SYSTEM, CONTINU-FLO SOLUTION SET, 0.2 MICRON HIGH PRESSURE EXTENDED LIFE FILTER, 2 INJECTION SITES, MALE LUER LOCK ADAPTER
2C8571¹	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET (10 DROPS/ML), 0.2 MICRON HIGH PRESSURE EXTENDED LIFE FILTER, 2 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C8593¹	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, 0.2 MICRON DOWNSTREAM FILTER, 3 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
ACT8570¹	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, 0.2 MICRON DOWNSTREAM FILTER, LUER ACTIVATED VALVE, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR

CODE	BRIEF DESCRIPTION
------	-------------------

SPECIALTY SETS

1C8043	NITROGLYCERIN SET WITH DUO-VENT SPIKE, PVC TUBING SEGMENT, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR
2C7551 ¹	INTERLINK SYSTEM, NITROGLYCERIN SET WITH DUO-VENT SPIKE, PVC TUBING SEGMENT, INJECTION SITE, MALE LUER LOCK ADAPTER
2C8851 ¹	CLEARLINK SYSTEM, NITROGLYCERIN SET WITH DUO-VENT SPIKE, PVC TUBING SEGMENT, LUER ACTIVATED VALVE, MALE LUER LOCK ADAPTER
2C7557 ^{1,2}	INTERLINK SYSTEM, VENTED PACLITAXEL SET (10 DROPS/ML), POLYETHYLENE LINED TUBING, NON-DEHP PUMP SEGMENT, 0.2 MICRON FILTER INJECTION SITE, MALE LUER LOCK ADAPTER
2C7558 ^{1,2}	INTERLINK SYSTEM, PACLITAXEL SET, POLYETHYLENE LINED TUBING, NON-DEHP PUMP SEGMENT, 0.2 MICRON FILTER INJECTION SITE, MALE LUER LOCK ADAPTER
2C8857 ^{1,2}	CLEARLINK SYSTEM, VENTED PACLITAXEL SET (10 DROPS/ML), POLYETHYLENE LINED TUBING, NON-DEHP PUMP SEGMENT, 0.2 MICRON FILTER, LUER ACTIVATED VALVE, MALE LUER LOCK ADAPTER
2E7554 ¹	SOLUTION SET, MALE EPIDURAL NON-LUER LOCK ADAPTER. FOR EPIDURAL USE. NRFIT [TM].

STOPCOCK AND MANIFOLD SETS

2C6254	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, 2 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR. LARGE BORE STOPCOCK MANIFOLD EXTENSION SET, LUER ACTIVATED VALVE, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR.
2C6255	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, 2 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER. LARGE BORE 4-WAY STOPCOCK EXTENSION SET, 2 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER.
2C6256	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, 2 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR. LARGE BORE 4-WAY STOPCOCK MANIFOLD EXTENSION SET, LUER ACTIVATED VALVE, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR.
2C6511	INTERLINK SYSTEM, CONTINU-FLO SOLUTION SET, 2 INJECTION SITES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR. LARGE BORE STOPCOCK MANIFOLD EXTENSION SET, INJECTION SITE, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR.
2C6931	INTERLINK SYSTEM, CONTINU-FLO SOLUTION SET, 2 INJECTION SITES, MALE LUER LOCK ADAPTER, 3-PORT MANIFOLD WITH CHECK VALVES. LARGE BORE 4-WAY STOPCOCK EXTENSION SET, 2 INJECTION SITES, MALE LUER LOCK ADAPTER.
2C8931	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, 2 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER. 3-PORT MANIFOLD WITH CHECK VALVES, 1.0 ML. LARGE BORE 4-WAY STOPCOCK EXTENSION SET WITH ROTATING MALE LUER LOCK, 2 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR.
3C0062	INTERLINK SYSTEM, CONTINU-FLO SOLUTION SET, 2 INJECTION SITES, MALE LUER LOCK ADAPTER. LARGE BORE 4-WAY STOPCOCK EXTENSION SET, 2 INJECTION SITES, MALE LUER LOCK ADAPTER.
3C0091	INTERLINK SYSTEM, CONTINU-FLO SOLUTION SET, 2 INJECTION SITES, MALE LUER LOCK ADAPTER. STOPCOCK EXTENSION SET, LARGE BORE STOPCOCK WITH ROTATING MALE LUER LOCK. EXTENSION SET, INJECTION SITE, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR.
3C0105	INTERLINK SYSTEM, CONTINU-FLO SOLUTION SET, 2 INJECTION SITES, MALE LUER LOCK ADAPTER. LARGE BORE STOPCOCK WITH ROTATING MALE LUER LOCK EXTENSION SET, INJECTION SITE, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR. EXTENSION SET, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR.
3C0134	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, 2 LUER LOCK ADAPTER WITH RETRACTABLE COLLAR. LARGE BORE 4-WAY STOPCOCK MANIFOLD EXTENSION SET, LUER ACTIVATED VALVE, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR.
3C0162	CLEARLINK SYSTEM, CONTINU-FLO SOLUTION SET, 2 LUER ACTIVATED VALVES, MALE LUER LOCK ADAPTER. LARGE BORE STOPCOCK MANIFOLD EXTENSION SET, LUER ACTIVATED VALVE, MALE LUER LOCK ADAPTER. EXTENSION SET, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR.
STP0073	INTERLINK SYSTEM, SOLUTION SET, INJECTION SITE, MALE LUER LOCK ADAPTER. LARGE BORE STOPCOCK WITH ROTATING MALE LUER LOCK EXTENSION SET, INJECTION SITE, MALE LUER LOCK ADAPTER WITH RETRACTABLE COLLAR.

¹Indicates an administration set is designated as an OTHER compatible set.

²Indicates a set with set change interval of 96 hours or 3 L, whichever comes first.

This compatible administration set list is subject to change over time.
Please visit <https://service.baxter.com/documents> for the most current version.

Baxter.com

Baxter, Clearlink, Continu-Flo, Duo-Vent, Interlink and Novum IQ are trademarks of Baxter International Inc. or its subsidiaries.

US-MD34-240012 v1.0 03/2024